

# Vital Sign Monitor




## BT-720

- 4.3" color touch screen
- SpO2, Pulse, NIBP, Masimo SpO2
- Short/long trend graphic / trend table display
- Audio & visual alarm, adjustable alarm limit
- Automatic brightness adjustment
- Accurate SpO2 performance during motion and low perfusion
- Perfusion Index data / bar graph display
- Pitch tone variation for pulse rate
- Internal memory for data storage
- SD card for easy software upgradation
- Over 8 hours continuous working on rechargeable lithium-ion battery
- Specialized PC software for data review and analysis
- Option for Masimo SpO2, NIBP
- HL7 export to clinical information systems



Portable design




SpO2




SpO2+NIBP



Trend Graphic



Trend Tabular



PC Viewer Software

# Pulse Oximeter




## BT-710

- 4.3" color touch screen
- Handheld style
- SpO2, Pulse
- Accurate SpO2 performance during motion and low perfusion
- Specialized PC software for data review and analysis
- Over 8 hours continuous working on rechargeable lithium-ion battery
- 5V DC input with convenient Micro-USB charger
- Battery can be charged by external USB battery
- Suitable for adult, pediatric and neonate



Handheld size



Trend Tabular



## BT-700

- 3.2" color LCD
- Desktop Pulse Oximeter
- SpO2, Pulse
- Trend Graphic, Trend Tabular display
- Internal data saving
- Firmware upgrade by USB interface
- Apply nurse call state